

Small Size Low-Speed Centrifuge TD4

Features

1. Small Size and light weight with great space saver for the lab.
2. The small LCD display, shows the Speed, RCF, Time Acc and Dec at same times.
3. High performance and low noise;
4. The base is fixed by sucker, suitable for vehicle.
5. Free maintenance DC brush-less motor drive, running stable and quiet.
6. Necessary equipment for the small hospital clinics, and laboratories..



Parameters

Model	TD4		
Max Speed	4000r/min		
Rotor	Capacity	Max Speed	Max RCF
No.1 Angle Rotor	8*15ml	4000r/min	1880xg
Motor	Microprocessor & Brushless DC Motor		
Timer Range	0~99min		
Noise	≤50dB(A)		
Power Supply	AC220V&110V 50Hz 2A		
Dimension	320*235*206mm(L*W*H)		
Weight	5kg		