

Benchtop Low-Speed Centrifuge TD4A-WS

TD4A-WS is suitable for routine sample analysis in medical, hospital, pathology and institutional laboratories. With wide variety of accessories, they can be also used for preparation of samples in industrial and research laboratories.

Features:

1. Adopt the steel structure, which is safer, simple-compact, light-weight and low-noise.
2. Digital display and LCD display for choose, adopt control system realizing microprocessor control, it control rotate speed, relative centrifuge force.
3. Maintenance-free frequency conversion motor, can prevent over-speed operation.



Parameters

Model	TD4A-WS		
Speed	0-4000rpm		
No.1 Angle Rotor	4/6×50ml	4000r/min	2680xg
No.2 Swing Rotor	4×50ml	4000r/min	2680xg
No.3 Swing Rotor	12×10ml	4000r/min	2680xg
No.4 Swing Rotor	8×10/15ml	4000r/min	2680xg
Motor	Brushless DC Motor, microprocessor		
Display	Digital display / LCD display		
Speed Accuracy	±10rpm		
Noise	≤50dB(A)		
Timer	0~99min		
Dimension	450*340*270mm		
Weight	19kg		