

Micro Hematocrit Centrifuge TG12-W

TG12-W capillary blood centrifuge is mainly used to determine the value of the blood hair cells, as well as the centrifuge of the micro blood and micro capacity.

Features

- **Stainless steel chamber**, compact structure.
- **Safety Interlock system** that prevents opening the centrifuge lid during running.
- **With Emergency lid lock release**, in case of power off.



Parameters

Model	TG12-W
Max speed	12000r/min
Max RCF	14800xg
Rotor Capacity	24 Pieces of Capillary
Timer Range	0-99min
Motor	Microprocessor control, DC Brush-less motor
Noise	≤55dB(A)
Power Supply	AC220V&110V 50Hz 2A
Dimension	475*355*245mm(L*W*H)
Weight	17.5KG