

Micro Low Speed Refrigerated Centrifuge TDL4-MC

Introduction

Micro Low Speed Refrigerated Centrifuge for small capacity tubes special.

Features

- Advanced CPU microprocessor control, brushless motor.
- Compressor temperature control and keep low temperature working environment.
- All parts of rotor can be High Temperature and High Pressure Sterilization.
- LCD display , emergency door lid open, and running stable with low noise , low vibration.
- With Angle Rotor, high efficiency and compact structure.



Parameters

Model	TDL4-MC		
NO.1 Angle Rotor	8*15ml	4000r/min	3500xg
Motor	Microprocessor & Brushless DC Motor		
Temperature Range	-20°C~+40°C		
Timer Range	0~99min		
Noise	≤50dB(A)		
Power Supply	AC220V&110V 50Hz 2A		
Dimension	540*340*312mm(L*W*H)		
Weight	40kg		