

## Tabletop PRP Centrifuge TD4-PRP

### Introduction

TD4-PRP is Fat and PRP purified centrifuge, it can use 10ml, 20ml syringe to separate directly, the safe transfer process ensure the high survive rate of PRP and fat ,and reduced the complex operation process, it is the idea machine for cosmetic surgery, beauty salon, clinics and so on . can also fit many type of PRP kits in market.

### Features

1. Achieved the maximum amount of PRP extraction through clinical research and test and special design on the Acc and Dec and RCF.
2. Microprocessor control, Frequency conversion brushless AC motor.
3. Touch LCD display with PRP program, easy to operate.
4. Protection from imbalance, lid-lock, over-speed and display the error code on the screen.
5. It is equipped with free time brake and able to set time, speed, Acc/Dec according to the requirement.



### Parameters

Model	TD4-PRP	
Max Speed/RCF	4000r/min / 3050xg	
NO.1 Swing Rotor	No.1 Swing Rotor (SY20)	4×20ml syringe (with bucket cap)
NO.2 Swing Rotor	No.2 Swing Rotor (PX60)	4×60/30ml PRO PRP tube
NO.3 Swing Rotor	No.3 Swing Rotor (PX100)	4×20ml ACP double syringe
NO.4 Swing Rotor	No.4 Swing Rotor (PX20)	4×20ml DR PRP Tube
NO.5 Swing Rotor	No.5 Swing Rotor (DS15)	4×15ml ACP double syringe (with bucket cap)
NO.6 Swing Rotor	No.6 Swing Rotor (DS20)	4×20ml ACP double syringe (with bucket cap)
Motor	Microprocessor & Brushless DC Motor	
Timer Range	0~99min	
Noise	≤50dB(A)	
Power Supply	AC220V&110V 50Hz 2A	
Dimension	470*395*320mm(L*W*H)	
Weight	29kg	