

Tabletop PRP & FAT Centrifuge TD4-ZF

Introduction

TD4-ZF is multi-functional PRP & self fat purification centrifuge, which is mainly used in cosmetic surgery, dermatology and TCM rehabilitation the urology field.

It can uses 10ml, 20ml, 30ml,50ml of conventional syringes to separate the fat & PRP and no heating or damage on the fat, easier to operate and ensured survive rate of the fat.



Parameters

Model	TD4-ZF		
Max Speed	5000r/min		
Max RCF	4390xg		
Speed Accuracy	±20r/min		
NO.1 Swing Rotor	4*50/20ml syringe	4000r/min	4390xg
No.2 Swing Rotor	4*30ml Tricell PRP Tube	4000r/min	3500xg
Motor	Microprocessor & AC Motor		
Timer Range	0~99min		
Noise	≤50dB(A)		
Power Supply	AC220V&110V 50Hz 5A		
Dimension	545*452*370mm(L*W*H)		
Weight	39kg		