

## Tabletop PRP/CGF/FAT Centrifuge TDL4-ZF

### Introduction

TDL4-ZF has 7-inch touch screen, built-in special preparation programs for PRP, PRF, and CGF, can also set parameters freely. Adopt fluorine-free high-efficiency refrigeration system(-20°C~+40°C). 10 kinds of lifting and lowering gears. Equipped with door cover protection, overspeed and electronic unbalance detection, protection system and UV lamp automatic sterilization device.



### Parameters

<b>Model</b>	<b>TDL4-ZF</b>	
<b>Max Speed</b>	<b>5000r/min</b>	
<b>Max RCF</b>	<b>5220xg</b>	
<b>Rotor</b>	<b>No.1 Swing Rotor (SY20)</b>	<b>4×20ml syringe (with bucket cap)</b>
	<b>No.2 Swing Rotor (PX60)</b>	<b>4×60/30ml PRO PRP tube</b>
	<b>No.3 Swing Rotor (PX100)</b>	<b>4×20ml ACP double syringe</b>
	<b>No.4 Swing Rotor (PX20)</b>	<b>4×20ml DR PRP Tube</b>
	<b>No.5 Swing Rotor (DS15)</b>	<b>4×15ml ACP double syringe (with bucket cap)</b>
	<b>No.6 Swing Rotor (DS20)</b>	<b>4×20ml ACP double syringe (with bucket cap)</b>
<b>Temp Accuracy</b>	<b>±1°C</b>	
<b>Temp Range</b>	<b>-20°C~+40°C</b>	
<b>Timing Range</b>	<b>1sec-99min59sec</b>	
<b>Motor</b>	<b>Microprocessor &amp; AC Motor</b>	
<b>Noise</b>	<b>≤60dB(A)</b>	
<b>Power Supply</b>	<b>AC220V 50Hz 5A</b>	
<b>Dimension</b>	<b>485*550*360mm(L*W*H)</b>	
<b>Weight</b>	<b>45kg</b>	