

DW-23 Monocular Student Microscope



Introduction

This series microscope is a basic student microscope. It has new design, reasonable price, optimum performance, novel construction and sharp image. It is suitable for teaching in primary and middle schools and also for human medical routine examination in hospitals.

Specifications

Wide Field Plan-scope Eyepiece	WF10X/18mm; Total Magnification: 40X-640X
Achromatic Objectives	4X, 10X, 40X(S)
Head	Monocular Head, 45° Inclined, fixed head
Single Stage with Clips	95X95mm
Single Lens Condenser	N.A.0.65; Six-stop Dial Diaphragm
Coarse Focusing Adjustable Mechanism	8mm
Illumination	LED Lamp
Optional	1.WF16X, WF20X; 2.Rechargeable Batteries 3.Top LED Lamp
Net Weight	1.5kg

Package

1set/Foam + Inner Carton
6sets/Outer Carton, Gross Weight: 12KG
Carton Dimension: 55×49×38cm