

## DW-23A Monocular Student Microscope



### Introduction

This series microscope is a basic student microscope. It has new design, reasonable price, optimum performance, novel construction and sharp image. It is suitable for teaching in primary and middle schools and also for human medical routine examination in hospitals.

### Specifications

Wide Field Plan-scope Eyepiece	WF10X/18mm;Total Magnification: 40X-640X
Achromatic Objectives	4X,10X,40X(S)
Head	Monocular Head, 45°Inclined,360° Rotatable
Double layers Mechanical Stage	95X105mm, moving range: 50X20mm
Single Lens Condenser	N.A.0.65; $\Phi$ 2- $\Phi$ 30mm Iris Diaphragm
Coaxial Coarse Focusing Adjustable Mechanism	17mm
Illumination	LED Lamp
Optional	1.WF16X, WF20X; 2.Rechargeable Batteries 3.Top LED Lamp
Net Weight	1.5kg

### Package

1set/Foam + Inner Carton

6sets/Outer Carton, Gross Weight: 12KG

Carton Dimension: 55×49×38cm