

DW-2005B Binocular Biological Microscope



Introduction

This series microscope is a high precision instrument designed and produced by our company. It is especially designed for clinical examination and teaching demonstration in medical and health establishments, laboratories, agricultural science and technology field, research institutes. It is used for routine work and research in biology, bacteriology cytology and pharmacology.

Specifications

Wide Field Plan-scope Eyepiece	WF10X
45mm DIN Achromatic Objective	4X,10X,40X(S), 100(S,Oil)
Articulated Free Binocular Head, Interpupillary Distance	55-75mm; 30°Inclined, 360°Rotatable
Double Layers Mechanical Stage	125X125mm, Moving Range X-Y: 60X40mm
ABBE Condenser	N.A.1.25; Φ 2- Φ 30mm Iris Diaphragm and Φ 32 Filter
Coaxial Coarse and Fine Focusing Adjustable Mechanism	14mm
Illumination	Built-in Adjustable Brightness 1W LED lamp
Optional	1.WF10X, WF16X; 2.Semi Plan Objective, Plan Objective; 3.Rechargeable Batteries
Net Weight	4kg

Package

1set/Foam + Inner Carton

2sets/Outer Carton;Gross Weight/Carton: 10.8kg

Carton Dimension: 44×37×47cm