

DW-2102B Binocular Biological Microscope



Introduction

This series microscope is a bigger, basic model and high precision instrument designed and produced by our company. It is especially designed for clinical examination and teaching demonstration in medical and health establishments, laboratories, agricultural science and technology field, research institutes. It is used for routine work and research in biology, bacteriology cytology and pharmacology

Specifications

Wide Field Plan-scope Eyepiece	WF10X/18mm
45mm DIN Achromatic Objective	4X,10X,40X(S), 100(S,Oil)
Head	Articulated Free Binocular Head, Interpupillary Distance:55-75mm;30° Inclined , 360°Rotatable
Double Layers Mechanical Stage	130X140mm, Moving Range X-Y: 70X50mm
ABBE Condenser	N.A.1.25; Φ 2- Φ 30mm Iris Diaphragm and Φ 32 Filter
Coaxial Coarse and Fine Focusing Adjustable Mechanism	20mm, Precision:0.002mm
Illumination	Built-in Adjustable Brightness 1W LED lamp; 6V/20W Halogen Lamp
Optional	1.WF10X, WF16X; 2.Semi Plan Objective, Plan Objective; 3.Rechargeable Batteries 4.Electronic Eyepiece 1.3Mp, 3.0Mp
Net Weight	5.6kg

Package

2sets/Foam + Carton

Carton Dimension: 53×37×49cm; Gross Weight: 14kg