

DW-9106E Monocular Student Microscope



Introduction

This series microscope is a basic model and high precision instrument designed and produced by our company. It is especially designed for clinical examination and teaching demonstration in medical and health establishments, laboratories, agricultural science and technology field, research institutes. It is used for routine work and research in biology, bacteriology cytology and pharmacology.

Specifications

Wide Field Plan-scope Eyepiece	WF10X/18mm; Total Magnification: 40X-1600X
45mm DIN Achromatic Objective	4X,10X,40X(S),100X(S,Oil)
Head	Monocular Head, 45°Inclined, 360°Rotatable
Double Layers Mechanical Stage	110x125mm, Moving Range X-Y:60x20mm
ABBE Condenser	N.A. 1.25; Φ 2- Φ 30mm Iris Diaphragm And Φ 32 Filter
Coaxial Coarse and Fine Focusing Adjustable Mechanism	14mm
Illumination	6V/20W Halogen Lamp, LED Lamp
Optional	1.WF16X, WF20X; 2.Semi Plan Objective, Plan Objective; 3.Rechargeable Batteries
Net Weight	2.6kg

Package

1set/Foam + Inner Carton
4sets/Outer Carton, gross weight: 13KG
Carton Dimension: 54×39×40cm