

# PO10 Portable pH/ORP Meter

## Features

- 3.7” Color screen, user-friendly interface.
- pH Calibration points: 1-3 points.
- ORP: measure and display temperature (°C/°F) simultaneously.
- Data logging up to 3000 groups.
- Equipped with a USB Type-C charging port, conforms with the EU new law.
- Ideal for use in harsh environments — IP57 Waterproof and Dustproof rating, silicone socket protector and grip cover protect the instrument from general impact.
- Large backlit LCD displays measuring value and temp along with icons of stable readings.
- Foldable stand allows use as a benchtop meter.
- For general liquid test.



## Specifications

pH	Range	(-2.000-20.000) pH
	Resolution	0.01/0.1/0.001 pH
	Accuracy	±0.002 pH ± one digit
ORP (mV)	Range	-1,999mV - +1,999mV
	Resolution	± 0.1/1mV
	Accuracy	±0.03% FS ±1digit
Temperature	Range	0-100°C
	Resolution/Accuracy	0.1°C; ±0.2°C±1digit
Data storage/content		3000 groups (Serial number, Mode, Measuring value, Temp, Date: Y/M/D/H/M)
Data export/Charging port		Type-C - USB
PC software function		GLP measurement management Two-level user management The user's operation history log Meter usage management Automatically /manually timing data storage
Power supply		Rechargeable Lithium battery × 1
Packing including		1. Meter× 1 2. ORP301 plastic ORP probe× 1 3. Temp-02H Stainless steel probe× 1 4. 222mV ORP standard solution, 50mL× 1 5. Plastic ATC pH electrode 201T-H× 1 6. pH buffer solution: 4.00/7.00/10.01 50mL, one for each 7. Type-C Charger Cable× 1 8. Hard carrying case× 1