

DW-SI-II Benchtop Shaking Incubator


Features

- ◆ Integrated with incubator and shaker to save space and cost.
- ◆ Polished stainless-steel chamber, powder coated steel housing.
- ◆ LCD screen to display temperature and shaking speed.
- ◆ With operation data memory function to eliminate faulty operation.
- ◆ Auto-recovery after power outage or crash as originally programmed, avoiding data loss.
- ◆ Auto-stop operation when door opened.
- ◆ Brushless DC motor, more stable and reliable.



Optional Accessories

1. Universal spring shaking platform
2. Flask Clamps



Specifications

Model	DW-SI-103B-II	DW-SI-100B-II	DW-SI-200B-II
Volume	26L	48L	69L
Shaking Speed	20-300 rpm		
Speed accuracy	±1 rpm		
Swing amplitude	Φ26 mm		
Temp. Range	RT+5~60℃	RT+5~60℃	4~60℃
Temp. Uniformity	±1℃		
Temp. Accuracy	±0.1℃		
Flask Clamps(standard)	100ml*9	50ml*4 100ml*4 250ml*3 500ml*3	50ml*5 100ml*5 250ml*4 500ml*3
Flask Clamps (optional)	50ml*12 or 100ml*9 or 250ml*9 or 500ml*4	50ml*20 or 100ml*9 or 250ml*12 or 500ml*9 or1000ml*4	100ml*20 or 250ml*16 or 500ml*9 or 1000ml*5 or 2000ml*4
Tray Size	295×253mm*1	400×370mm*1	450×410mm*1
Chamber Size (W×D×H)	370*320*220mm	440×405×270mm	480*450*320mm
Exterior Size(W×D×H)	560*480*440mm	685*625*610mm	830*720*640mm
Timing Range	1-9999min		