

DW-LVO Vacuum Drying Oven

Features

- ★ Microprocessor controller with LCD display, more accurate and reliable.
- ★ Polished stainless-steel chamber, durable and easy to clean.
- ★ Designed with a transparent window made of reinforced tempered glass, easy observation.
- ★ Door tightness can be adjusted, silicon sealing.
- ★ Equipped with leakage protection and a spare mechanical thermostat to ensure the product work normally even the main temperature control failed.

Options

1. Programable temperature controller
2. Built-in Printer
3. RS485 Connector
4. Inert gas connection



Specifications

Model	DW-LVO-6020	DW-LVO-6050
Volume	25L	54L
Temp. Range	RT+10~200℃	RT+10~250℃
Temp. Fluctuation	±1℃	
Display Resolution	±0.1℃	
Vacuum Degree	<133 Pa	
Power Supply	AC200V, 50Hz	
Power Rating	0.5KW	1.4KW
Chamber Size (W×D×H)	30×30×28cm	42×35×37cm
Exterior Size (W×D×H)	58×57×45cm	71×53×59cm
Solid Racks (Std.)	1	2
Heating Technology	Shelfing heating	Jacket heating
Vacuum Pump	Optional	Optional