

# DW-LVO Floor Standing Vacuum Drying Oven

## Features

- ★ Microprocessor controller with LCD display, more accurate and reliable.
- ★ Polished stainless-steel chamber, durable and easy to clean.
- ★ Designed with a transparent window made of reinforced tempered glass, easy observation.
- ★ Independent temperature control for each solid shelf.
- ★ Door tightness can be adjusted, silicon sealing.
- ★ Equipped with leakage protection.

## Options

1. Programable temperature controller
2. Built-in Printer
3. RS485 Connector
4. Inert gas connection



## Specifications

Model	DW-LVO-6090	DW-LVO-6210	DW-LVO-6500	DW-LVO-6933
Volume	91L	215L	428L	1007L
Temp. Range	RT+10~200℃			
Temp. Fluctuation	±1℃			
Display Resolution	±0.1℃			
Vacuum Degree	<133 Pa			
Power Supply	AC200V, 50Hz			
Power Rating	1.6KW	2.2KW	3.8KW	5.5KW
Chamber Size (W×D×H)	45×45×45cm	56×60×64cm	63×81×84.5cm	75.5×116×115cm
Exterior Size (W×D×H)	61×66×145cm	72×82×175cm	100×104×185.5cm	93×130×182cm
Heating Technology	Shelf heating			
Solid Racks (individual control)	2	3	4	5