

DW-LVO-LC

Vacuum Drying Oven with Pressure Control

Features

- ★ Microprocessor controller with pressure and temperature display.
- ★ Polished stainless-steel chamber, durable and easy to clean.
- ★ Designed with a transparent window made of reinforced tempered glass, easy observation.
- ★ Door tightness can be adjusted, silicon sealing.
- ★ Independent temperature control for each solid shelf.
- ★ Equipped with leakage protection.

Options

1. Multi-segment programmable control
2. Built-in printer
3. RS485 interface
4. Inert gas connection



Specifications

Model	DW-LVO-6050LC	DW-LVO-6090 LC	DW-LVO-6210 LC
Volume	54L	91L	215L
Temperature Range	RT+10~250 °C	RT+10~200°C	
Temperature Fluctuation	±1°C		
Display Resolution	0.1°C		
Vacuum Degree	<133 Pa		
Power Supply	AC 220V, 50Hz		
Power Rating	1.4 KW	1.6 KW	2.2 KW
Chamber Size (W×D×H) cm	42×35×37cm	45×45×45cm	56×60×64cm
Exterior Size (W×D×H)	71×53×137cm	61×66×145cm	72×82×175cm
Heating Technology	Jacket heating	Shelf heating	Shelf heating
Solid Shelves	2	2(individual ctrl)	3(individual ctrl)