

# DW-DWT Water bath

## Features

- ◆ PID control, more accurate.
- ◆ Built-in temperature calibration and LCD display, easy operation.
- ◆ Seamless stainless-steel chamber and powder coated steel housing. Safe and easy to clean.
- ◆ With power leakage protection.
- ◆ With timing function, timing range: 1-9999 minutes.



## Specifications

Model	DW-DWT-600	DW-DWT-420
Temperature Range	RT+5°C to 95°C	
Temperature Range	±0.5°C	
Temperature Uniformity	±0.5°C	
Power Supply	220V／50Hz	
Power Rating （W）	1800±10%	900±10%
Chamber Size (L*W*H) mm	505*300 *150	300*240*150
Volume	23L	11L
Exterior Size （L*W*H）mm	655*335*300	450*275*300
Net/Gross Weight	11.5/21kg	7.5/14kg