

Precision Balance 1mg

Features

- 1. Electromagnetic Loadcell
- 2. Ultra High Definition LCD
- 3. Multi Unit Conversion
- 4. RS232 Bidirectional Interface
- 5. Internal Calibration
- 6. Quick Sensitivity Setting
- 7. Anti-Static Glass Windshield



Specifications

External Cal.	DT123	DT223	DT323	DT523	DT1003	DT1003N
Internal Cal	DT123N	DT223N	DT323N	DT523N		
Accuracy class	Ⓜ	Ⓜ	Ⓜ	Ⓜ	Ⓜ	Ⓜ
Max Capacity	100g	200g	300g	500g	600g	1000g
Readability	0.001g					
Repeatability	±0.001g					
Linearity	±0.002g					
Stable time	≤3s					
Operation Temperature	15℃~30℃, Fluctuation ≤ 2℃/h					
Relative Humidity	40%~80%					
Units	g/ct/oz/mg					
Scale Pan Size	Φ 90mm					
Product Size	320*210*320mm					
Windshield Size	225*190*265mm					
Gross Weight	7.2kg					
Power	Input 220V 50Hz Output DC9V-2A					