

FA Ecnomic Balance 0.1mg

Features

- 1. Electromagnetic Loadcell
- 2. Ultra High Definition LCD
- 3. Multi Unit Conversion
- 4. RS232 Bidirectional Interface
- 5. Internal Calibration
- 6. Quick Sensitivity Setting
- 7. Anti-Static Glass Windshield



Specifications

External Cal.	FA1204	FA2204	FA3204
Internal Cal	FA1204N	FA2204N	FA3204N
Accuracy class	Ⓢ	Ⓢ	Ⓢ
Max Capacity	120g	220g	320g
Readability	0.0001g		
Repeatability	±0.0001g		
Linearity	±0.0002g		
Stable time	≤3s		
Operation Temperature	15℃～30℃, Fluctuation ≤ 2℃/h		
Relative Humidity	40%~80%		
Units	g/ct/oz/mg		
Scale Pan Size	Φ 90mm		
Product Size	320*210*320mm		
Windshield Size	225*190*265mm		
Gross Weight	7.2kg		
Power	Input 220V 50Hz Output DC9V-2A		