

DW-CMII Clarity Tester

Description

DW-CMII type clarity tester adopts visual light inspection method, which is used for the clarity test of all kinds of injections, I.V. fluids and bottled medicinal liquids.

Features

- ◆ The instrument uses three-colour fluorescent lamps, electronic ballasts and adjustable shading devices composed of the light path system, which can eliminate stroboscopic, adjustable illumination, improve the resolution of visual inspection, reduce visual fatigue.
- ◆ Digital electronic illuminance meter reading intuitive, stable and reliable. Detection time can be set arbitrarily with sound and light alarm function.
- ◆ Detection time can be set arbitrarily, and has sound and light prompt function. Tips timing start and end.
- ◆ The instrument is reasonable in design, stable in performance, easy in operation and beautiful in appearance, and is the ideal instrument for clarity testing as stipulated in the Chinese pharmacopoeia.



Specifications

Power supply	AC 220V±10%, 50Hz
Power(Output)	22W
Illumination Range	1000~6000lx
Work time	1~99S
Size	703mm×175mm×513mm
Weight	9.4Kg