

DW-CSD01 Digital Friability Test Apparatus

Description

Tablets are inevitably subjected to vibration or friction during production, transportation, etc. These factors may cause breakage of tablets and affect application. Tablet friability is an indicator of a tablet's ability to resist shock and abrasion.It is generally measured using a tablet friability tester.



Specifications

Time setting	0~4minutes
Rotation count	100 turns
Timing Range	0--9:59:00
Counting Range	0--99999
RPM	25±1times/minutes
Power(output)	20W
Power supply	AC 220V 50HZ
Dimensions	440*286*350mm
Dial size	insidedia286mm, depth39mm
Net weight	10KG