

SHZ-III Circulating Water Vacuum Pump

Features

The pump is on the basis of circulating water multi-purpose vacuum pump, according to the characteristics of the small laboratory area, improve the transparent anti-corrosion tank, anti-corrosion shell and stainless steel shell, reduce the size of the improvement, with small volume, light weight, beautiful appearance and other characteristics, double table, double suction head.

Features

- The machine adopts two taps and can be used separately or in parallel with two vacuum gauges.
- The main machine adopts two types of stainless steel movement and anti-corrosion material movement.
- Corrosion resistance, no pollution, low noise, easy to move, can also be installed according to user needs vacuum regulating valve.



Specifications

Model	SHZ-III
Power	180W
Working power supply	220V/50Hz
Flow rate	60L/min
Head	10M
Material	Anti-corrosion engineering plastics
Maximum vacuum	-0.098Mpa
Single head pumping capacity	10L/min
Number of pumping heads	2
Tank capacity	15L
External dimensions	400*280*420mm
Pump pressure	0.3bar
Weight	10.5kg