

SHZ-95B Circulating Water Vacuum Pump

Features

- 5 exhaustors can be used individually or in parallel, with a vacuum gauge.
 - Pump body is made of stainless steel; motor made by Sino-foreign joint venture; pump spindle made of stainless steel; Pump head design complying with fluid mechanics; shell is made of engineering plastics or stainless steel.
- The exhaust speed is as fast as 14L/MIN.
- Anti-corrosive, without pollution, low noise, easy to move.
 - Optional: vacuum regulating valve.
 - This vacuum pump can absorb the harmful gas, which is environmental friendly and ensures the physical and psychological health of experimenters.



Specifications

Model	SHZ-95B
Power	550
Voltage	220V/50Hz
Flow Rate	≥80L/min
Lift	≥12m
Safety Protection	Check valve old rolled sheet with antistatic spray
Shell material	Electrostatic spray cold plate
Max. Vacuum Degrr	-0.098MPa
Unit Extraction	10L/min
Exhaustor	5
Water tank volume	≥50L
Water tank material	PVC
Weight	35kg