

TF Series **Electronic Balance 0.01g**



Introductions

High precision strain gauge loadcell, skip the tedious operation process. Dual interface , dual display optional, provide more intelligent applications in the field of strain gauges, and provide customers with good cost performance and convenient use experience.

Features

- High contrast LCD display
- Glass windshield
- Unit conversion, counting, percentage, etc.
- RS232C interface (Option USB interface)
- Dual display (Optional)
- Big size stainless steel scale pan
- Automatic fault detection, overload protection
- Tare function
- Below weighing facility (Optional)
- 4 x AA battery (Optional)

Specifications

Model	TF2002 TF2002S	TF3002 TF3002S	TF5002 TF5002S	TF6002 TF6002S	TF10002 TF10002S	TF20002 TF20002S	TF30002 TF30002S	TF50002 TF50002S
Capacity	0-200g	0-300g	0-500g	0-600g	0-1000g	0-2000g	0-3000g	0-5000g
Readability	0.01g							
Repeatability(≤)	±0.02g							
Linearity(≤)	±0.02g							
Interface	RS232C/USB(Optional)							
Pan Size	170×170mm							
Product Size	32×19.5×7cm							
Packing Size	24.4×12.1×39.5cm							
N.W.	1200g							
G.W.	1500g							
Power	110V/220-240V AC							