

TF Series Electronic Balance 0.1g



Introductions

High precision strain gauge loadcell, skip the tedious operation process. Dual interface , dual display optional, provide more intelligent applications in the field of strain gauges, and provide customers with good cost performance and convenient use experience.

Features

- High contrast LCD display
- Glass windshield
- Unit conversion, counting, percentage, etc.
- RS232C interface (Option USB interface)
- Dual display (Optional)
- Big size stainless steel scale pan
- Automatic fault detection, overload protection
- Tare function
- Below weighing facility (Optional)
- 4 x AA battery (Optional)

Specifications

Model	TF2001 TF2001S	TF6001 TF6001S	TF10001 TF10001S	TF15001 TF15001S	TF20001 TF20001S	TF30001 TF30001S	TF50001 TF50001S
Capacity	0-200g	0-600g	0-1000g	0-1500g	0-2000g	0-3000g	0-5000g
Readability	0.1g						
Repeatability(≤)	±0.2g						
Linearity(≤)	±0.2g						
Interface	RS232C/USB(Optional)						
Pan Size	170×170mm						
Product Size	32×19.5×7cm						
Packing Size	24.4×12.1×39.5cm						
N.W.	1200g						
G.W.	1500g						
Power	110V/220-240V AC						