

DW-MW Series Halogen Moisture Analyzer



Features

1. Halogen lamp/HBM sensor
2. Aluminum shell/Stainless steel chamber
3. LCD(White back light with black font)
4. Temperature/Time set
5. The percentage of moisture content/dry residual
6. Store historical set
7. Option: Aluminum box/Printer/Grinder/USB



Specifications

Model	DW-110MW	DW-105MW	DW-102MW	DW-100MW
Capacity	110g			
Readability	0.01g	0.005g	0.002g	0.001g
Operating Temp	5°C-35°C			
Display	LCD(White back light with black font)			
Calibration	External Calibration			
Interface	RS232 (standard) / USB (optional)			
Pan Size	Φ90mm			
Packing Size	490x350x360mm			
Gross Weight	7.5kg			
Heating Source	Halogen Lamp-1			
Temp.Sensor	PT100			
Temp.Range	40°C-199°C			
Temp,Step	1°C			
Moisture Range	0.09%-100.0%	0.009%-100.00%		
Moisture Readability	0.2%	0.10%	0.04%	0.01%
Dry Rang	100.00%-0.00%			
Dry Readability	0.2%	0.10%	0.04%	0.01%
Temp Setting	40°C-199°C By 1°C Step			
Time Setting	1-99min By 10s			
Number of Storage	15			
Heater	220V±15% 50HZ/110V±15% 60HZ/400W			
Balance	Input: 220V 50HZ/110V 60HZ Output:9V			
Chamber Height	22mm			