

DW-F800-B Cold Extractor



Characteristics

- DW-F800-B can simultaneously de-fat 6 samples before hot extraction, solvent dehydration of fiber residues and lignin determination.
- DW-F800-B use the same crucible system as DW-F800, permitting samples to be dried and weighed between extractions if required.
- Equipped with high-pressure pump to prevent sample plug crucible, guarantee the extraction efficient.
- With a splash guard, eliminate potential safety problems of highly corrosive liquid during the extraction.

Specifications

Working Temperature	10-40°C
Sample Weight	0.5g-3g
Storage Temperature	10-50°C
Batch Capacity	6 pcs per batch
Relative Humidity	0-90%
Power Supply	220 VAC $\pm$ 10% 50Hz